

METHOD OF USE

INSTRUCTION FOR USE

Distributed by / Distribuidor:
Coloplast Corp
1601 West River Road North
Minneapolis, MN 55411 USA
www.coloplast.us
Tel: 1-800-533-0464

Australian sponsor / Distributed by
Coloplast Pty Ltd
Level 4, 1 Acacia Place,
Fernree Business Park,
Notting Hill, VIC. 3168
Australia
www.coloplast.com.au
Freecall Numbers
Australia 1800 653 317
New Zealand 0800 600 018

Distributed by / Distribué par:
Coloplast Canada Corporation
1380 Creditstone Road, Unit 6&7
Concord, ON L4K 0J1
Tel: 1-888-880-8605

Rev. 04-2019

Plastod S.p.A.
40012 Lippo di Calderara
Bologna, Italy
www.plastod.com

Distributor:
Coloplast A/S
Holtevej 1
3050 Humlebaek,
Denmark

ENGLISH

Silicone contact layer for atraumatic removal and wound protection

PRODUCT DESCRIPTION

Biatain Contact is a sterile biocompatible silicone contact layer, for atraumatic removal and wound protection.

Easy to apply, it stays in place and must be used in combination with a secondary absorbent dressing, which needs to be fixated.

It is composed of a perforated transparent polyurethane support film and of a wound direct contact layer in silicone gel. Its biocompatible composition, together with the highly conformable structure and open holes, allows the drainage of moisture and exudates from the wound to the secondary dressing, ensures high breathability and permeability to exudates, reducing the risk of stagnation and tissue maceration.

The structure limits the risk of cell ingrowth within the pores and on the dressing surface, while favoring the tissue granulation process of the wound bed. Its characteristics of conformability, together with the high transparency allow to monitor the healing process for the whole duration of the application.

It protects the wound bed even at secondary dressing changes, as it prevents adhesion to newly formed tissues. It minimises possible trauma for both the wound and perilesional skin, limiting the risk of pain on removal. It prevents the adhesion of the secondary dressing to the wound, leaving no residue on removal.

INDICATIONS OF USE

Biatain Contact is a contact dressing, indicated for the treatment of ulcers and chronic or acute lesions, at various levels of exudation, to be treated with a concomitant application of a secondary absorbent dressing.

It can also be used as a protective layer, after appropriate fixation with adequate dressing, on non-exuding wounds or on fragile skin.

On specific medical indications, it can also be used with compressive bandages.

METHOD OF USE

Biatain Contact offers the advantage of staying in place during the application and to be easily removed at the end of use.

Please follow the procedure below:

Clean, disinfect and medicate lesion and skin according to the protocol of reference.

Choose Biatain Contact in the size sufficient to cover the entire lesion and the perilesional skin. If necessary, Biatain Contact can be cut to the most suitable size using sterile scissors.

- 1 Remove the first protective liner, exposing the silicone layer.
- 2 Gently apply the silicone layer of Biatain Contact, cover completely the area of the lesion, whilst removing the second protective liner.
- 3 Apply a sterile and fastening absorbent dressing on Biatain Contact.

CHANGE FREQUENCY

The frequency of dressing change depends on the status of the wound and in accordance with the current good clinical practice. Monitor the wound. In case of complications or infections, seek for medical advice.

CONTRAINDICATIONS

Do not use Biatain Contact in patients with known sensitivity to silicone.

PRECAUTIONS AND SYMBOLS

Do not reuse (single use). Reuse may cause alteration in performance and cross-contamination.

Sterile. Sterilized by Ethylene Oxide.

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